


MICROCNX® SERIES CONNECTORS

The **MicroCNX® Standard and Luer** are micro-connectors that provide a simple, efficient method of connecting tubing for small-format biomanufacturing assemblies. The MicroCNX series is engineered specifically for the challenging conditions of biological media transfer in bioprocessing, cell therapy, and gene therapy applications. The connector comes in a range of tubing sizes featuring a luer connection option to create a sterile, closed connection on an existing luer connection.



SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar (Standard)
Up to 75 psi, 5.1 bar for 48 hours (Standard)
Up to 43.5 psi, 3 bar (Luer)

OPERATING TEMPERATURE

34°F to 104°F (1°C to 40°C)

STORAGE TEMPERATURE

-112° F to 140° F (-80° C to 60° C)

TERMINATIONS

1/16", 3/32", 1/8" ID hose barb (1.6mm, 2.4mm, 3.2mm), Male Luer; Female Luer

MATERIALS

Main Components: Polycarbonate (white)

Seals: Silicone (clear), platinum-cured

Protective Cover: Polypropylene (teal),

Membrane: Hydrophobic Polyethersulfone

STERILIZATION

Gamma: Up to 50kGy irradiation.

Autoclave: One cycle up to 266°F (130°C) for 60 minutes

FEATURES

PINCH-CLICK-PULL →

Easy-to-Use →

Luer adapter →

BENEFITS

Easy to use three-step connection process reduces risk of operator error

No hassling with external machines or equipment, or mismatched tubing

Connect to a luer/open format connection in a biosafety cabinet, remove from the biosafety cabinet, and have a sterile connection anywhere in your production.

TYPICAL FLOW RATE

Cv Value Range: 0.04-0.27

Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

DID YOU KNOW

MicroCNX connectors eliminate the need to purchase, calibrate, validate, maintain, and allocate clean room space for tube welding equipment.

Scan code to visit webpage



cpcworldwide.com/MicroCNX



COLDER PRODUCTS COMPANY

U.S.A.

PHONE: +1 (651) 645-0091

TOLL FREE: +1 (800) 519-7633

E-MAIL: info@cpcworldwide.com

COLDER PRODUCTS COMPANY GMBH

Germany

E-MAIL: cpcgmbh@cpcworldwide.com

DOVER (SHANGHAI) INDUSTRIAL CO., LTD

Shanghai, China

PHONE: +86 21 2411 2666

TOLL FREE: +86 400 990 1978

E-MAIL: asiapacific@cpcworldwide.com

MICROCNX® SERIES DIMENSIONS

POLYCARBONATE with teal cover For autoclave or gamma irradiation applications.




MicroCNX® Standard



MicroCNX® Luer

TERMINATION	METRIC EQ.	PART NO.	A	B	C
1/16" HOSE BARB	1.6 mm	CNX17001HT	0.93" (23.6 mm)	2.05" (52.1 mm)	1.49" (37.8 mm)
3/32" HOSE BARB	2.4 mm	CNX17003HT	0.93" (23.6 mm)	2.16" (54.9 mm)	1.49" (37.8 mm)
1/8" HOSE BARB	3.2 mm	CNX17002HT	0.93" (23.6 mm)	2.25" (57.2 mm)	1.49" (37.8 mm)
FEMALE LUER		CNX170LFHT	0.93" (23.6 mm)	2.15" (54.6 mm)	1.49" (37.8 mm)
MALE LUER		CNX170LMHT	0.93" (23.6 mm)	2.15" (54.6 mm)	1.49" (37.8 mm)

PRODUCT DIMENSIONS



MICROCNX® CONNECTORS ASSEMBLY PROCEDURE

PINCH



Remove the protective cover on each half.

CLICK



Join halves and click together.

PULL



Pull membrane strips directly away from the connector.

WARRANTY: All sales are subject to Colder Products Company's limited express warranty set forth in the CPC catalog. Contact your local distributor or CPC Customer Service for warranty provisions.

Warning: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC. It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer.

COPYRIGHT © 2025 BY COLDER PRODUCTS COMPANY.

CPC, Colder Products Company, and Colder Products are registered trademarks with the United States Patent and Trademark Office. For detailed trademark information, please visit: <https://www.cpcworldwide.com/Trademarks>