

CONNECTION SOLUTIONS FOR BIOPHARMACEUTICAL PROCESSES

CPC A DOVER[®] COMPANY

WE INSPIRE CONFIDENCE AT EVERY POINT OF CONNECTION

CPC (Colder Products Company) is the leader in the design and manufacture of single-use connection technology and connectors for the biopharmaceutical market. CPC offers a wide variety of solutions including sterile connect, sterile disconnect, SIP connections and open connects. Our innovative designs provide flexibility for biopharmaceutical manufacturers to easily combine multiple components, single-use or hybrid systems including process containers, tubing manifolds, transfer lines, bioreactors, and other bioprocessing equipment.

Easy to use and robust single-use connectors from CPC maintain flow path sterility and integrity while enabling biopharmaceutical manufacturers to improve production yields, decrease time to market and reduce costs. Our genderless sterile connectors simplify process integration, maximize flexibility, and streamline supply chains. CPC makes peoples' lives better by developing innovative high-quality products that make media transfer safe and easy.

QUALITY

At CPC, everyone is involved in meeting or exceeding our customers' expectations from our suppliers to our distribution network, and most important, our employees. CPC measures and continually improves our standards of product quality, support services and overall customer and employee satisfaction. CPC's Quality System conforms to ISO 9001 and ISO 13485 standards. Products for biopharmaceutical applications are manufactured in our two ISO Class 7 certified cleanrooms.

Learn more about our cleanrooms and quality control processes at: cpcworldwide.com/resources-support/quality-compliance



EXPERIENCE

CPC is the leading provider of quick connect couplings, fittings, disconnects and combination connectors used in fluid transfer. We innovate, engineer and manufacture fit-for-purpose products in close collaboration with our customers.

Founded in Minnesota in 1978, CPC has built a successful, growing company by focusing exclusively on critical points of connection within fluid management systems. CPC's reach is global with operations in the US, Germany and China, sales offices in ten countries and hundreds of distributor partners and OEM solution providers around the world.

CPC's biopharma team includes an innovative research and development group solely focused on creating the next generation of single-use technologies, with dedicated product managers, technical specialists, quality engineers and test lab expertise.

Our single-use and closed systems connectors empower our customers' solutions to be safer, more efficient and reliable. CPC biopharma experts provide media handling expertise for our customers as well as to the industry. These resources include:

- Channel Management Team supporting our OEM and integrator partners.
- Applications Development Team serving as consultants to our end user customers.
- Customer Fulfillment and Inside Sales Teams serving our channels and end customers.

CPC offers application and operator training designed to provide guidance on where and how single-use technology can be used or optimized in the manufacturing process. For more details visit cpcworldwide.com/training.

SUPPLY CHAIN

With a commercial model designed to be market neutral, CPC is committed to offering consistent pricing programs, common delivery lead times and product availability information to all customers. Our market neutrality ensures that CPC retains supply chain integrity, ensures accountability and upholds our reputation as a world-class organization.

Additionally, as part of CPC's commitment to meeting the needs of biopharma customers, our products are produced in multiple cleanroom manufacturing facilities. This redundancy is designed to maintain product availability, manufacturing efficiency and reliability of manufacturing processes.

ONLINE RESOURCES

Visit cpcworldwide.com/bio for answers to your questions about our company and products.

VALIDATION TEST REPORTS

Validation test reports provide details of all the testing that has been performed on the product to ensure confidence at every point of connection. Extractables data can also be requested.

PRODUCT VIDEOS

Check out some of the latest innovations, technologies, and product tutorials in our CPC biopharma videos. Our video library contains instructions on how to assemble CPC connectors as well as best practices and tips to ensure you get the most out of every connection.

APPLICATION ARTICLES

CPC's industry experts share their knowledge on specific biopharmaceutical applications and how single-use technologies have helped our customers improve production.

CAD MODELS

See our connectors from all angles, anytime, anywhere. Get immediate access to 2D drawings and 3D models for use in your manufacturing process diagrams or to fit your specific needs. Simply register on our site to download CAD models in many different file formats.

REGULATORY & COMPLIANCE DOCUMENTS

CPC follows strict regulatory compliance standards to ensure the quality of our supply of our products. The materials used in our broad portfolio of products are compliant with various regulatory bodies including NSF, RoHs, REACH and more. Download the documents you require to meet regulatory and compliance standards.

ASK OUR ENGINEERS

We're here to help. When you have questions, CPC's team of expert engineers has answers. From flow rates to material compatibility and more, we specialize in providing media fluid connection solutions to meet the requirements of your most complex biopharma applications. Looking for something more specific to your needs, reach out to one of our industry experts with a specific question.

VISIT CPCWORLDWIDE.COM/BIO



REGULATORY AND COMPLIANCE

ISO 13485 CERTIFICATION

ISO 13485 is recognized by regulators around the world as a good basis for addressing medical device design and manufacturing regulatory requirements. It allows us to enhance product safety by proactively identifying and managing product and project risks. Our quality management system is ISO 13485 certified, which allows us to better control the consistency of manufactured products.

ISO 9001 CERTIFICATION

ISO 9001 is a standard which assures consistency of a product ordered by customers. Organizations having ISO 9001 certification have demonstrated compliance to the ISO 9001:2015 requirements by an independent registration authority. CPC's Quality Management System has been approved and certified under the ISO 9001 standard.

CLEANROOM MANUFACTURING

CPC manufactures certain Life Sciences and Chemical Management product lines in cleanrooms certified by an external testing service to meet or exceed ISO Class 7 at 0.5 mm per ISO 14644. Certification data is available upon request.

ANIMAL DERIVED COMPONENT FREE (ADCF)

According to declarations from CPC's raw material suppliers, the materials used to manufacture the flow path components of the biopharmaceutical product lines do not contain substances of animal origin.

FDA

The U.S. Food and Drug Administration publishes, through the Code of Federal Regulations, standardized criteria which govern the acceptability of materials used in food contact.

REACH

REACH is the regulation for Registration, Evaluation, Authorisation and Restriction of Chemicals. It entered into force on 1st June 2007 to streamline and improve the former legislative framework on chemicals of the European Union (EU). REACH places greater responsibility on industry to manage the risks that chemicals may pose to the health and the environment. CPC publishes a list of CPC products that are compliant with the EU regulation 1907/2006.

RESTRICTION OF HAZARDOUS SUBSTANCES (ROHS)

This directive bans new electrical and electronic equipment containing more than agreed levels of lead, cadmium, mercury, hexavalent chromium, polybrominated biphenyl (PBB) and polybrominated diphenyl ether (PBDE) flame retardants.





UNDERSTANDING SINGLE-USE SYSTEMS

Increasing global demand for new biologics, vaccines and cell therapies is driving manufacturers to replace traditional stainless equipment with single-use systems, which consist of plastic-based processing equipment used in the development and production of biopharmaceutical drugs.

BENEFITS OF SINGLE-USE

Operational Efficiencies	Increases flexibility and faster batch turnaround.							
Cost Effectiveness	Minimizes cleaning and validation requirements.							
Economic Advantages	Reduces capital expenditures and facility footprints.							
Safety and Quality	Improves sterility assurance while decreasing the risk of cross-contamination and product loss.							
Flexibility	Facilitates multi-drug production and fast product changeover.							
Sustainability	Consumes less water, energy and chemicals when compared to stainless-based processing. Single-use plastic waste is an excellent fuel source for waste-to-energy conversion.							



UNDERSTANDING THE BIOPHARMACEUTICAL MANUFACTURING PROCESS







Single-use bioprocessing is designed to be flexible, efficient, and effective in the manufacture of drug substances, monoclonal antibodies, vaccines, biosimilars and regenerative medicines. The process is split into two main sections, upstream and downstream.

Along with growing the cell line, the purpose of the upstream process is to scale-up the volume of the target protein or cell from volumes as small as a vial to bioreactors that can be as large as 5,000L.

Once the target specimen has reached a target yield, it turns the corner towards the downstream process. The goal of the downstream line is to clarify, purify and filter the target. Reaching this goal is done with processes such as clarification, viral inactivation, chromatography and various types of filtration. The final step is fill and finish, where the target cell has been grown and purified to the point where it can be used for filling syringes or other devices to give to a patient to start their healing process.

ASEPTIC CONNECTION TECHNOLOGY* PAGES 10-23

12 MICROCNX® SERIES CONNECTORS: Sterile connection for your 1/8" flow applications

MATERIAL: Polycarbonate, silicone TERMINATIONS: 1/16"; 3/32" & 1/8" ID hose barb (1.6mm, 2.4mm and 3.2mm)

16 ASEPTIQUIK® S SERIES CONNECTORS: Sterile connection for your 1/4" flow applications MATERIAL: Polycarbonate, silicone

TERMINATIONS: 1/8"; 1/4"; 3/8" ID hose barb (3.2mm, 6.4mm and 9.5mm); 1/4"; 3/4" sanitary & MPC insert

18 ASEPTIQUIK[®] G SERIES CONNECTORS: Sterile connection for your 1/2" flow applications MATERIAL: Polyphenylsulfone, polycarbonate and silicone TERMINATIONS: 1/4"; 3/8"; 1/2"; 3/4" ID hose barb (6.4mm, 9.5mm, 12.7mm, 19.0mm) & 3/4"; 1 1/2" sanitary & MPC insert

20 ASEPTIQUIK[®] L SERIES CONNECTORS: Sterile connection for your 1" flow applications MATERIAL: Polycarbonate, silicone TERMINATIONS: 3/4"; 1" ID hose barb (19.0 mm; 25.4 mm); 1-1/2" sanitary

22 ASEPTIQUIK® W SERIES CONNECTORS: Sterile connection for your 1 1/2" flow applications

MATERIAL: Polycarbonate, silicone TERMINATIONS: 1"; 1 ¼"; 1 ½" ID hose barb (25.4mm; 31.75mm; 38.1mm); 1-1/2" sanitary

STERILE DISCONNECTION TECHNOLOGY* PAGES 24-29

28 HFC DISCONNECT SERIES

MATERIAL: Polysulfone, alloy spring, and silicone TERMINATIONS: 1/4"; 3/8" & 1/2" ID hose barb (6.4mm, 9.5mm & 12.7mm)

ASEPTIC COMBINATION CONNECTORS PAGES 30-35

34 ASEPTIQUIK[®] G DC SERIES CONNECTORS: Sterile connection and disconnection combination for 1/2" flow applications

MATERIAL: Polycarbonate

TERMINATIONS: 1/4"; 3/8"; 1/2" ID hose barb (6.4mm; 9.5mm and 12.7mm)















OPEN FORMAT CONNECTION TECHNOLOGY PAGES 36-49

40 MPC SERIES CONNECTORS: Open-format connection for your 3/8" flow applications MATERIAL: Polycarbonate, polysulfone and silicone TERMINATIONS: 1/8"; 1/4" & 3/8" ID hose barb (3.2mm, 6.4mm, 9.5mm)

42 MPX SERIES CONNECTORS: Open-format connection for your 1/2" flow applications MATERIAL: Polycarbonate, polysulfone, and silicone TERMINATIONS: 3/8"; 1/2" ID hose barb (9.5mm & 12.7mm)

MPC/MPX BACK-TO-BACK SERIES ADAPTERS: Connect single-use systems that may feature identical connections at the end of their tubing. MATERIAL: Polycarbonate, polysulfone and silicone

46 MPC/MPX SANITARY SERIES CONNECTORS: Attaches directly to 3/4", 1" & 1-1/2" sanitary terminations to provide greater flexibility for integrating components into single-use or hybrid process systems.

MATERIAL: Polysulfone and silicone TERMINATIONS: 3/4", 1"; 1-1/2" sanitary

44

48 MPU SERIES CONNECTORS: Open-format twist-to-lock connection for 3/4" flow applications

MATERIAL: Polysulfone and silicone TERMINATIONS: 3/4" (19.0mm) & 1" ID (25.4 mm)

STEAM-IN-PLACE CONNECTION TECHNOLOGY PAGES 50-57

54 STEAM-THRU[®] SERIES CONNECTORS: Hybrid connection between stainless steel and single use applications

MATERIAL: Polysulfone and silicone

TERMINATIONS: 3/8" & 1/2" ID hose barb (9.5mm & 12.7mm) 3/4"; 1 1/2" sanitary

56 ASEPTIQUIK[®] STC SERIES CONNECTORS: Hybrid connection between stainless steel and single use applications, with a single-use AseptiQuik included.

MATERIAL: Polycarbonate, polysulfone and silicone TERMINATIONS: 3/4"; 1 1/2" sanitary

*If you don't see a product line, please contact your CPC representative for more information.

















ASEPTIC CONNECTION TECHNOLOGY

MICROCNX® SERIES CONNECTORS

MicroCNX® Connectors introduce a new category of aseptic micro-connectors that provide a simple, efficient method of connecting tubing for small-format biomanufacturing assemblies. MicroCNX connectors are the modern alternative to the cumbersome, industrial process of tube welding. Building on the inventiveness of CPC, the leader in single-use connection technology, the MicroCNX line of connectors is engineered specifically for the challenging conditions of biologic media transfer in bioprocessing, cell therapy and gene therapy applications.



SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours

OPERATING TEMPERATURE 34°F to 104°F (1°C to 40°C)

TERMINATIONS 1/16", 3/32" 1/8" ID hose barb (1.6mm, 2.4mm, 3.2mm)

MATERIALS

Main Components: Polycarbonate (white) Seals: Silicone (clear), platinum-cured Protective Cover: Polypropylene (teal) Membrane: Hydrophobic Polyethersulfone

STERILIZATION

Gamma: Up to 50kGy irradiation. **Autoclave:** One cycle up to 266°F (130°C) for 60 minutes

FEATURES

CPC Click -

TYPICAL FLOW RATE

Cv Value Range: 0.04-0.27

for MicroCNX hose barb terminations



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Caution: CPC does not recommend the use of PVC tubing with the polycarbonate MicroCNX[®] Series Connectors. Leachables from PVC tubing could be incompatible with polycarbonate material and affect product performance. Determining product application suitability is solely the customer's responsibility. CPC does not guaranty or warrant product suitability for any application or use.

BENEFITS

Intuitive three-step connection process reduces risk of operator error

Lowers risk of operator error and related performance, reliability and safety concerns

Eases single-use systems specifications with one part number for both halves

Audible confirmation of assembly with no additional hardware required

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/MicroCNX

DID YOU KNOW

MicroCNX connectors eliminate the need to purchase, calibrate, validate, maintain, and allocate clean room space for tube welding equipment.

Scan code to visit webpage



cpcworldwide.com/MicroCNX

MICROCNX® SERIES DIMENSIONS

POLYCARBONATE with teal cover *For autoclave or gamma irradiation applications*.



TERMINATION 1/16" Hose Barb	METRIC EQ. 1.6 mm	PART NO. CNX17001HT	▲ 0.93" (23.6 mm)	B 2.05" (52.1 mm)	C 1.49" (37.8 mm)
3/32" HOSE BARB	2.4 mm	CNX17003HT	0.93" (23.6 mm)	2.16" (54.9 mm)	1.49" (37.8 mm)
1/8" HOSE BARB	3.2 mm	CNX17002HT	0.93" (23.6 mm)	2.25" (57.2 mm)	1.49" (37.8 mm)

PRODUCT DIMENSIONS



MICROCNX® CONNECTORS ASSEMBLY PROCEDURE



Caution: CPC does not recommend the use of PVC tubing with the polycarbonate MicroCNX[®] Series Connectors. Leachables from PVC tubing could be incompatible with polycarbonate material and affect product performance. Determining product application suitability is solely the customer's responsibility. CPC does not guaranty or warrant product suitability for any application or use.



WHAT IS ASEPTIQUIK® STERILE TECHNOLOGY?

AseptiQuik® Connectors provide quick and easy sterile connections, even in non-sterile environments–a critical capability for biopharmaceutical manufacturers. Featuring a straightforward, simple three-step connection process and a wide range of termination options—including 1/8- to 1 1/2-inch sizes and genderless connections—the AseptiQuik series allows you to transfer media easily with less risk of error. Their robust, reliable performance eliminates the need for clamps, fixtures or tube welders, giving you sterile, high-quality single-use connections every time.



FEATURES

BENEFITS

Genderless design	\longrightarrow	Eases integratio universal mating series
Rubust construction	\longrightarrow	Repeatable and additional hardv
FLIP-CLICK-PULL	\longrightarrow	Innovative three- of operator error
Integrated pull tab covers	\longrightarrow	Pull tabs act as complexity and e membranes
CPC Click	\longrightarrow	Audible confirma

ases integration of single-use systems with niversal mating between connectors of the same eries

Repeatable and reliable performance with no additional hardware required

Innovative three-step connection process reduces risk of operator error (AQS, AQG, & AQL)

Pull tabs act as protective cover reducing part complexity and ensure simultaneous removal of both membranes

Audible confirmation of assembly

NOMINAL FLOW PATH (NFP) SIZE



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the NFP.

ASEPTIQUIK CONNECTORS OPTIONS

	TERMINATIONS								MATERIALS		STERILIZATION (CHOOSE 1)		FLOW INFORMATION						
	HOSE BARB					SANITARY Flanges			Le la	-BOMAL	ALL SHIT	OMPHIBIL	Hell is	HHHHHHHHHHHHHH					
	1,0%	"h	200	1st	214	"/	1/1911	I'll	1 th	10 m	I'mer	MD SHI		A COL	AL POIL	St. Shimp		A AONINAL	C WALLE
ASEPTIQUIK® S SERIES	\checkmark	\checkmark	\checkmark						\checkmark	\checkmark		\checkmark	\checkmark	\checkmark		\checkmark	\checkmark	1/4"	0.19-1.74
ASEPTIQUIK® G SERIES		~	\checkmark	\checkmark	\checkmark					\checkmark	\checkmark	\checkmark	\checkmark	\checkmark		\checkmark	\checkmark	1/2"	1.5 - 31
ASEPTIQUIK® G PPSU SERIES		\checkmark	\checkmark	\checkmark	\checkmark					\checkmark			~		\checkmark	\checkmark		1/2"	1.5 - 31
ASEPTIQUIK® L SERIES					\checkmark	\checkmark					\checkmark		\checkmark	\checkmark		\checkmark	\checkmark	1"	30-57
ASEPTIQUIK® W Series						\checkmark	\checkmark	\checkmark			\checkmark		~	\checkmark		~	\checkmark	1 1/2"	39-118

ASEPTIQUIK CONNECTORS ASSEMBLY PROCEDURE



PULL

To complete the connection, simply snap the pull tabs together by pushing on the CPC logos and pull the membranes from the AseptiQuik connector halves.



CLICK



Scan code to watch AseptiQuick assembly video



https://youtu.be/un2PnvUAZ0w

Note – For AQW assembly see product page

ASEPTIC STERILE CONNECTION

ASEPTIQUIK[®] S Series connectors

AseptiQuik® S Connectors provide quick and easy sterile connections for small-flow applications, even in non-sterile environments. The "FLIP-CLICK-PULL" design of AseptiQuik S enables users to easily transfer small volumes of media with less risk of operator error than with traditional methods. The connector's genderless and robust design provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers can now make 1/8", 1/4" and 3/8" hose barb and 1/4" and 3/4" sanitary sterile connections with the quality and market availability they expect from the leader in single-use connection technology.

SPECIFICATIONS

OPERATING PRESSURE Up to 60 psi, 4.1 bar

OPERATING TEMPERATURE 39°F to 104°F (4°C to 40°C)

STERILIZATION

Gamma: Up to 50kGy irradiation AutoClave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

TERMINATIONS

 $1/8",\,1/4"$ and 3/8" ID hose barb (3.2mm, 6.4mm and 9.5mm), 1/4" and 3/4" sanitary and MPC insert

MATERIALS

Main Components: Polycarbonate (white) Pull Tabs/Caps: Polycarbonate (blue, standard version) Polycarbonate (white, HT version) Seals: Silicone (clear), platinum-cured Membrane:

Polyethylene (standard version) Hydrophobic polyethersulfone (HT version), PTFE strip sticker

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES



AQS ¼" Sani with Smooth Bore

TYPICAL FLOW RATE:

Cv Value Range: 0.19 - 1.74 for AseptiQuik S



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/AseptiQuik-S



BENEFITS

Eases single-use systems specifications with onepart number for both halves

Intuitive three-step connection process reduces risk of operator error

Audible confirmation of connection with no additional hardware required

Enables the ability to change a BPC or other singleuse system with open format connections to closed systems

Minimizes transitional flow disruptions throughout upstream processing

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-S

DID YOU KNOW

Did you know that the AseptiQuik S is perfect for simplifying the process to pull flask samples from your bioreactor?



ASEPTIQUIK S SERIES DIMENSIONS

POLYCARBONATE with blue pull tabs - For gamma irradiation applications.



TERMINATION	METRIC EQ.	PART NO.		B	G
1/8" HUSE BAKB	3.2 mm	AQ517002	2.25° (57.2 mm)	1.30° (33.0 mm)	0.50° (12.7 mm)
1/4" HOSE BARB	6.4 mm	AQS17004	2.25" (57.2 mm)	1.45" (36.8 mm)	0.65" (16.5 mm)
1/4" SANITARY	6.4 mm	AQS33004	2.25" (57.2 mm)	1.50" (38.1 mm)	0.70" (17.8 mm)
3/8" HOSE BARB	9.5 mm	AQS17006	2.25" (57.2 mm)	1.45" (36.8 mm)	0.65" (16.5 mm)
3/4" SANITARY	19.1 mm	AQS33012	2.25" (57.2 mm)	1.60" (40.6 mm)	0.80" (20.3 mm)
MPC INSERT		AQS17MPC	2.25" (57.2 mm)	1.49" (37.9 mm)	0.69" (17.5 mm)

POLYCARBONATE HT with white pull tabs - For autoclave or gamma irradiation applications.

A	
	2
	8

TERMINATION	METRIC EQ.	PART NO.	А	B	C
1/8" HOSE BARB	3.2 mm	AQS17002HT	2.25" (57.2 mm)	1.30" (33.0 mm)	0.50" (12.7 mm)
1/4" HOSE BARB	6.4 mm	AQS17004HT	2.25" (57.2 mm)	1.45" (36.8 mm)	0.65" (16.5 mm)
1/4" SANITARY	6.4 mm	AQS33004HT	2.25" (57.2 mm)	1.50" (38.1 mm)	0.70" (17.8 mm)
3/8" HOSE BARB	9.5 mm	AQS17006HT	2.25" (57.2 mm)	1.45" (36.8 mm)	0.65" (16.5 mm)
3/4" SANITARY	19.1 mm	AQS33012HT	2.25" (57.2 mm)	1.60" (40.6 mm)	0.80" (20.3 mm)
MPC INSERT		AQS17MPCHT	2.25" (57.2 mm)	1.49" (37.9 mm)	0.69" (17.5 mm)

PRODUCT DIMENSIONS



NOTES

ASEPTIQUIK® G Series connectors

AseptiQuik® G Connectors enable quick and easy sterile connections, even in non-sterile environments. The easy-to-use genderless design simplifies system integration and minimizes the risk of operator error. The connectors' robust construction provides enhanced user confidence and reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers benefit from a full range of interchangeable 1/4" to 1-1/2" termination solutions with the quality and market availability they expect from the leader in single-use connection technology.

SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours (excluding AQGMPC)

OPERATING TEMPERATURE 34°F to 104°F (1°C to 40°C)

STERILIZATION

Standard (blue) and PPSU (purple) Version: Gamma: up to 50kGy

High Temperature (white) Version: Gamma: up to 50kGy

Autoclave: One cycle up to 266°F (130°C) for 60 minutes

TERMINATIONS

1/4", 3/8", 1/2", 3/4" ID hose barb (6.4mm, 9.5mm, 12.7mm, 19.0mm) and 3/4", 1-1/2" sanitary

MATERIALS

Main Components:

Polycarbonate (white), (standard and HT versions) Polyphenylsulfone (off white) (PPSU version) **Pull Tabs/Caps:**

Polycarbonate (blue, standard), (white, HT version), (purple, PPSU version)

Seals:

Silicone (clear), platinum-cured

Membranes:

Polyethylene (standard and PPSU versions), Hydrophobic polyethersulfone (HT version), PTFE strip sticker

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES



BPA-free



BENEFITS

Eases single-use systems specifications with one-part number for both halves

Intuitive three-step connection process reduces risk of operator error

Audible confirmation of connection with no additional hardware required

AseptiQuik PPSU enables genderless connection for a greater range of chemical applications, offering versatile connections across downstream processes with a pH range from 2 to 12

AseptiQuik PPSU meets a broader range of Application Requirements

TYPICAL FLOW RATE:

Cv Value Range: 1.5 - 31 for AseptiQuik G



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/AseptiQuik-G

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-G

DID YOU KNOW

Did you know that the AseptiQuik G is perfect for connecting different buffers to your chromatography skid?



ASEPTIQUIK G SERIES DIMENSIONS

POLYCARBONATE with blue pull tabs - For gamma irradiation applications.



TERMINATION	METRIC EQ.	PART NO.	A	в	С
1/4" HOSE BARB	6.4 mm	AQG17004	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
3/8" HOSE BARB	9.5 mm	AQG17006	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
1/2" HOSE BARB	12.7 mm	AQG17008	2.62" (66.6 mm)	1.95" (49.5 mm)	1.15" (29.2 mm)
3/4" HOSE BARB	19.1 mm	AQG17012	2.62" (66.6 mm)	2.36" (59.9 mm)	1.56" (39.6 mm)
3/4" SANITARY	19.1 mm	AQG33012	2.62" (66.6 mm)	1.66" (42.2 mm)	0.86" (21.8 mm)
1-1/2" SANITARY	38.1 mm	AQG33024	2.62" (66.6 mm)	2.86" (72.6 mm)	2.02" (51.3 mm)
MPC INSERT		AQG17MPC	2.62" (66.6 mm)	1.50" (38.1 mm)	0.70" (17.8 mm)

POLYCARBONATE HT with white pull tabs - For autoclave or gamma irradiation applications.



TERMINATION	METRIC EQ.	PART NO.	A	B	0
1/4" HOSE BARB	6.5 mm	AQG17004HT	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
3/8" HOSE BARB	9.5 mm	AQG17006HT	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
1/2" HOSE BARB	12.7 mm	AQG17008HT	2.62" (66.6 mm)	1.95" (49.5 mm)	1.15" (29.2 mm)
3/4" HOSE BARB	19.1 mm	AQG17012HT	2.62" (66.6 mm)	2.36" (59.9 mm)	1.56" (39.6 mm)
3/4" SANITARY	19.1 mm	AQG33012HT	2.62" (66.6 mm)	1.66" (42.2 mm)	0.86" (21.8 mm)
1-1/2" SANITARY	38.1 mm	AQG33024HT	2.62" (66.6 mm)	2.86" (72.6 mm)	2.02" (51.3 mm)
MPC INSERT		AQG17MPCHT	2.62" (66.6 mm)	1.50" (38.1 mm)	0.70" (17.8 mm)

POLYPHENYLSULFONE with purple pull tabs - *For gamma irradiation applications.*

TERMINATION	METRIC EQ.	PART NO.	А	B	C
1/4" HOSE BARB	6.4 mm	AQG17104	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
3/8" HOSE BARB	9.5 mm	AQG17106	2.62" (66.6 mm)	1.62" (41.2 mm)	0.82" (20.8 mm)
1/2" HOSE BARB	12.7 mm	AQG17108	2.62" (66.6 mm)	1.95" (49.5 mm)	1.15" (29.2 mm)
3/4" HOSE BARB	19.1 mm	AQG17112	2.62" (66.6 mm)	2.36" (59.9 mm)	1.56" (39.6 mm)
3/4" SANITARY	19.1 mm	AQG33112	2.62" (66.6 mm)	1.66" (42.2 mm)	0.86" (21.8 mm)

PRODUCT DIMENSIONS



ASEPTIQUIK® L Series connectors

AseptiQuik® L Connectors enable quick and easy sterile connections, in large-volume, high-flow production environments. The large-format, 3/4", 1" hose barb and 1-1/2" sanitary genderless design simplifies system integration and minimizes the risk of operator error. The connectors' robust construction provides reliable performance without the need for clamps, fixtures or tube welders. Biopharmaceutical manufacturers benefit from the interchangeable connection sanitary flow solutions for full-scale bioprocessing production environments with the quality and market availability they expect from the leader in single-use connection technology.

SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours

OPERATING TEMPERATURE 34°F to 104°F (1°C to 40°C)

STERILIZATION

Gamma: Up to 50kGy irradiation Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

TERMINATIONS

3/4", 1" ID hose barb (19.0 mm, 25.4 mm) and 1-1/2" sanitary

MATERIALS

Main Components: Polycarbonate (white) Pull Tabs/Caps: Polycarbonate (blue, standard version) Polycarbonate (white, HT version) Seals: Silicone (clear), platinum-cured Membranes: Polyethylene (standard version) Hydrophobic polyethersulfone (HT version),

PTFE strip sticker

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES



FLIP-CLICK-PULL -----

CPC Click

Large Internal Diameter

TYPICAL FLOW RATE:

Cv Value Range: 30 - 57 for AseptiQuik L



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/AseptiQuik-L

BENEFITS

Eases single-use systems specifications with onepart number for both halves

Intuitive three-step connection process reduces risk of operator error

Audible confirmation of connection with no additional hardware required

Fast and efficient fluid transfer of large volumes under low pressures

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-L

DID YOU KNOW -

Did you know that the AseptiQuik L is perfect for connecting TFF, TFDF, and ATF and other filtration processes that require large flow volumes?

ASEPTIQUIK L SERIES DIMENSIONS

Polycarbonate with blue pull tabs - For gamma irradiation applications.



TERMINATION	METRIC EQ.	PART NO.	A	B	С
3/4" HOSE BARB	19.1 mm	AQL17012	4.22" (107.2 mm)	2.95" (74.9 mm)	1.50" (38.1 mm)
1" HOSE BARB	25.4 mm	AQL17016	4.22" (107.2 mm)	3.15" (80.0 mm)	1.70" (43.2 mm)
1-1/2" SANITARY	38.1 mm	AQL33024	4.22" (107.2 mm)	3.40" (86.4 mm)	1.95" (49.5 mm)

Polycarbonate HT with white pull tabs - For autoclave or gamma irradiation applications.

	TERMINATION	METRIC EQ.	PART NO.	A	B	C
~	3/4" HOSE BARB	19.1 mm	AQL17012HT	4.22" (107.2 mm)	2.95" (74.9 mm)	1.50" (38.1 mm)
4	1" HOSE BARB	25.4 mm	AQL17016HT	4.22" (107.2 mm)	3.15" (80.0 mm)	1.70" (43.2 mm)
1	1-1/2" SANITARY	38.1 mm	AQL33024HT	4.22" (107.2 mm)	3.40" (86.4 mm)	1.95" (49.5 mm)

PRODUCT DIMENSIONS





https://www.cpcworldwide.com/Blog/post/how-to-choose-sterile-connectors

ASEPTIC STERILE CONNECTION

ASEPTIQUIK[®] W SERIES CONNECTORS

Genderless AseptiQuik® W Series Connectors with nominal 1 1/2" flow path enable quick and easy sterile connections, in large-volume, high-flow production and process intensification environments. The genderless design simplifies system integration and minimizes the risk of operator error. The connectors' robust construction provides reliable performance without the need for clamps or fixtures. Biopharmaceutical manufacturers benefit from interchangeable 1", 1 ¼", and 1 ½" termination sizes for fullscale bioprocessing production environments with the quality and market availability they expect from the leader in the single-use connection technology.

SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours

OPERATING TEMPERATURE 34°F to 104°F (1°C to 40°C)

STERILIZATION

Gamma: up to 50 kGy irradiation Autoclave High Temp (HT) Version: Up to 266°F (130°C) for 60 minutes

TERMINATIONS

1", 1 1/4", 1 1/2" Hose barb, 1 1/2" Sanitary (25.4mm, 31.75mm, 38.1mm, 38.1mm)

MATERIALS

Main components: Polycarbonate (white) **Shipping Clip:** Polycarbonate (clear) Side latch: Polyphenylsulfone (blue) **Pull Tabs/Caps:** Polycarbonate (white, HT version) **Pull Tab Guide:** Polyphenylsulfone (gray) Seals: Silicone (clear), platinum-cured Membrane: Hydrophobic Polyethersulfone (HT version),

PTFE strip sticker

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES



BENEFITS

Eases single-use systems specifications with one-part number for both halves

Intuitive four-step connection process reduces risk of operator error

Increases connector robustness for wide variety of applications

Enables simultaneous and secure removal of both membranes

No need to purchase and store extra materials or tools to complete your sterile connection process

Move liquid medias faster to increase process efficiency

TYPICAL FLOW RATE:

Cv Value Range: 39-118 for AseptiQuik W



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-W



cpcworldwide.com/AseptiQuik-W

ASEPTIQUIK® W SERIES DIMENSIONS

POLYCARBONATE HT WITH WHITE PULL TABS - For autoclave or gamma irradiation applications.



TERMINATION 1" HOSE BARB	METRIC EQ. 25.4 mm	PART NO. Aqw17016ht	A 5.77" (146.6 mm)	B 4.00" (101.6 mm)	© 1.51" (38.4 mm)
1 1/4" HOSE BARB	31.8 mm	AQW17020HT	5.77" (146.6 mm)	4.39" (111.5 mm)	1.91" (48.5 mm)
1 1/2" HOSE BARB	38.1 mm	AQW17024HT	5.77" (146.6 mm)	4.91" (124.7 mm)	2.42" (61.5 mm)
1 1/2" SANITARY	38.1 mm	AQW33024HT	5.77" (146.6 mm)	3.66" (93.0 mm)	1.17" (29.7 mm)

PRODUCT DIMENSIONS



ASEPTIQUIK® W ASSEMBLY PROCEDURE





STERILE DISCONNECTION TECHNOLOGY

WHAT IS STERILE DISCONNECTION TECHNOLOGY?

CPC's sterile disconnection technology enables simple, quick sterile disconnection of your single-use systems with just the press of a thumb latch. Internal valves within the sterile disconnect technology close and seal upon disconnect, protecting the closed system fluid pathway on each side of the disconnected system. The HFC Disconnect product line enables sterile disconnections from tubing ¼" to ½" ID and connected sets feature a protective thumb latch to prevent accidental disconnections.



FEATURES

BENEFITS

Intuitive one-step disconnection process \longrightarrow	No requirement for additional equipment to make sterile disconnection Minimize operator error and ease standard operating procedure creation and training
Automatic shutoff values \longrightarrow	Stop flow upon disconnect
Protective thumb latch cover \longrightarrow	Eliminate accidental disconnects
Laser marked item number and lot number \longrightarrow	Full traceability to raw material source
Alloy C-276 internal flow path spring \longrightarrow	Enable broader application compatibility

OUR KNOWLEDGE IS YOURS.

Get useful information on SUT validation.



https://www.cpcworldwide.com/Biopharma-Campaigns/SUT-Validation-E-Book

NOMINAL FLOW PATH (NFP) SIZE



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the NFP.



HFC DISCONNECT OPTIONS

	TERMINATIONS		MATE	RIALS	PRODUCT Feature	STERIL (CHOC	IZATION DSE 1)	INFO	FLOW)RMATION	
	ŀ	IOSE BAR	B	Company and Com			MH COM	AND IN THE TOP TO THE T	Milling Cr.	Har wall
		3	11	S	23	8 1 1 B	Ch.	40	Ŵ	3
CONNECTED SET	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	3/8"	0.3 – 2.5
UNCONNECTED HALVES	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark		\checkmark	\checkmark	3/8"	0.3 – 2.5



HFC DISCONNECT SERIES CONNECTORS

HFC Disconnects enable sterile disconnection of singleuse biopharma and cell and gene therapy manufacturing systems. With an easy push of the connector thumb latch, sterility is maintained on both sides of the system during the disconnection process. The HFC Disconnect sets include protective thumb latch covers to help reduce the chance of accidental disconnection, and are laser marked with item and lot number for complete batch traceability.



SPECIFICATIONS

OPERATING PRESSURE Up to 75 psi, 5.17 bar

OPERATING TEMPERATURE 34° F to 104° F (1° C to 40°C)

STERILIZATION

Gamma: Up to 50 kGy irradiation - Autoclave: Up to 270°F (132°C), 60 minutes, one cycle

TERMINATIONS

1/4", 3/8" and 1/2" ID hose barb (6.4mm, 9.5mm and 12.7mm)

MATERIALS

Main components: Polysulfone (amber tint) **O-rings:** Silicone (clear), platinum-cured Flow Path Springs: Alloy C-276

WARNING: Pressure temperature chemicals and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

Scan code to visit webpage



cpcworldwide.com/HFC-Disconnect

FEATURES

Intuitive one-step disconnection process
Automatic shutoff valves
Protective thumb latch cover

Laser etched item number and lot number — Full traceability to raw material source

Alloy C-276 internal flow path spring —

\longrightarrow	No requirement for additional equipment to make sterile disconnection
\longrightarrow	Stop flow and eliminate need for pinch clamps
\longrightarrow	Guard against accidental disconnects

BENEFITS

 \longrightarrow Enabling broader application compatibility

TYPICAL FLOW RATE:

Cv Value Range: 0.3 - 2.5 for HFC Disconnect



Cy values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/HFC-Disconnect

DID YOU KNOW

The HFC Disconnect is great for post-use filter integrity testing (e.g. bubble point testing).

HFC DISCONNECT SERIES DIMENSIONS

COUPLING SETS - Polysulfone

6

	TERMINATION	METRIC EQ.	PART NO.	A	B
	1/4" HOSE BARB	6.4 mm ID	HFCD39SET4HC	1.54" (39.1 mm)	3.71" (94.3 mm)
AF EP	3/8" HOSE BARB	9.5 mm ID	HFCD39SET6HC	1.54" (39.1 mm)	3.71" (94.3 mm)
AN EN	1/2" HOSE BARB	12.5 mm ID	HFCD39SET8HC	1.54" (39.1 mm)	4.29" (109 mm)

COUPLING BODIES - Polysulfone

	TERMINATION	METRIC EQ.	PART NO.	Δ	в
C 100	1/4" HOSE BARB	6.4 mm ID	HFCD17439MHC	1.44" (36.6 mm)	2.78" (70.6 mm)
	3/8" HOSE BARB	9.5 mm ID	HFCD17639MHC	1.44" (36.6 mm)	2.78" (70.6 mm)
	1/2" HOSE BARB	12.5 mm ID	HFCD17839MHC	1.44" (36.6 mm)	3.07" (78.0 mm)

COUPLING INSERTS - Polysulfone

	TERMINATION	METRIC EQ.	PART NO.	A	в
Sin P	1/4" HOSE BARB	6.4 mm ID	HFCD22439MHC	1.00" (25.4 mm)	2.00" (50.8 mm)
COL N COL	3/8" HOSE BARB	9.5 mm ID	HFCD22639MHC	1.00" (25.4 mm)	2.00" (50.8 mm)
-	1/2" HOSE BARB	12.5 mm ID	HFCD22839MHC	1.00" (25.4 mm)	2.29" (58.2 mm)

MATING PARTS



PART	PART NO.	Α	8
SEALING CAP	HFC32039	1.44" (36.6 mm)	2.73" (69.3 mm)
SEALING PLUG	HFC30039M	1.00" (25.4 mm)	1.81" (46.0 mm)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters. Couplings are pictured with valves unless otherwise noted.

PRODUCT DIMENSIONS







WHAT IS ASEPTIC COMBINATION TECHNOLOGY?

Aseptic Combination Technology offers both a sterile connection AND a sterile disconnection in the same connector. This connect/disconnect combination connector helps process engineers and manufacturers make quick, easy, and sterile connections and disconnections, even in non-sterile environments. CPC's sterile combination technology is engineered using CPC's single-use technology, enabling users to make easy media transfers with less risk of operator error. When media transfers are complete, combination connectors easily and securely disconnect to help maintain sterility. Automatic shutoff valves support this disconnect, providing reliable performance without the need for tube welders.



FEATURES

n.		B. II.	-	-	-	~
ĸ.	the state	M	ы.	ы.		∕.
D		11				• •
_						-

${\sf Genderless} \longrightarrow$	Connects with any par product line
$FLIP \operatorname{-CLICK} \operatorname{-PULL} \longrightarrow$	Intuitive three-step co operator error
Simple Two-Step Disconnection \longrightarrow	Maintains system ster via disconnect process
Membrane Pull Tabs \longrightarrow	Ensure simultaneous a membranes
Automatic Shutoff Valves $\begin{tabular}{ c c c } \hline & & \\ \hline \hline & & \\ \hline & & \\ \hline \hline \\ \hline \hline & & \\ \hline \hline \hline \\ \hline \hline \\ \hline \hline \hline \\ \hline \hline \hline \\ \hline \hline \hline \hline \\ \hline \hline$	Eliminate the need for process
Protective Thumb Latch Cover \longrightarrow	Eliminates accidental
Alloy C-276 internal flow path spring \longrightarrow	Enables broader appli

connects with any part number within the AseptiQuik G roduct line ntuitive three-step connection process reduces risk of perator error Maintains system sterility in both sides of the system ia disconnect process insure simultaneous and secure removal of both nembranes liminate the need for tube sealers post disconnect rocess liminates accidental disconnects inables broader application compatibility

NOMINAL FLOW PATH (NFP) SIZE



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the NFP.



ASEPTIQUIK COMBINATION CONNECTORS OPTIONS

	TERMINATIONS								MA	TERIA	LS	STERIL (Choi	IZATION DSE 1)	I Info	FLOW RMATION	
			HOSE	BARB				SANITA Flang	RY ES	5	ht of	ALL AND	SUPPORT	Mahalt	41	ophilit and the second
	1,0%	1/4"	381	112	14	"/	1/411	341	I'm I'm	SILLOW	POLICE	2017PHIL	5 AMAR	Willogh's	HOMMAN	AND COMPANY
ASEPTIQUIK® G DC Series		\checkmark	\checkmark	\checkmark						~	\checkmark		\checkmark		1/2"	0.8 - 1.9

ASEPTIQUIK G DC ASSEMBLY PROCEDURE





Align the connector halves with the pull tabs hanging down. Slide the two halves together. Independently squeeze each side of the connector until you hear an audible "CPC Click."



To complete the connection, simply snap the pull tabs together by pushing on the CPC logos and pull the membranes from the connector.

ASEPTIQUIK G DC DISASSEMBLY PROCEDURE



CPC

ASEPTIQUIK® G DC SERIES CONNECTORS

AseptiQuik® GDC Series Connector is the first all-in-one, genderless, single-use connection technology to offer both a sterile connection and a sterile disconnection. With the AseptiQuik G DC connector, manufacturers can make quick and easy sterile connections and disconnections – even in non-sterile environments.

The intuitive Flip - Click - Pull design of the AseptiQuik G DC connector enables users to transfer media easily with less risk of operator error. After transfer is complete, a simple two-step disconnection maintains media sterility by preventing environmental ingress into the media flow path.

The connector's robust design and automatic shutoff valves provide reliable performance without the need for fixtures or tube welders/sealers.



SPECIFICATIONS

OPERATING PRESSURE

Up to 60 psi, 4.1 bar Up to 75 psi, 5.1 bar for 48 hours

OPERATING TEMPERATURE 34°F to 104°F (1°C to 40°C)

STERILIZATION Gamma: up to 50 kGy irradiation

TERMINATIONS

1/4", 3/8", 1/2" hose barb (6.4mm, 9.5mm, 12.7mm)

MATERIALS

Main components: Polycarbonate (white) Thumb latch: Polycarbonate (blue) Pull Tabs/Caps: Polycarbonate (blue) Seals: Silicone (clear), platinum-cured Membrane: Polyethylene Springs: Alloy C-276

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES

Genderless	\longrightarrow
FLIP - CLICK - PULL	\longrightarrow
Simple Two-Step Disconnection -	\longrightarrow
Membrane Pull Tabs	\longrightarrow
Automatic Shutoff Valves	\longrightarrow
Protective Thumb Latch Cover —	\longrightarrow
Alloy C-276 internal flow path spri	$_{1g} \longrightarrow$

TYPICAL FLOW RATE:

Cv Value Range: 0.8-1.90 for AseptiQuik G DC hose barb terminations



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

BENEFITS

operator error

- Connects with any part number within the AseptiQuik G product line Intuitive three-step connection process reduces risk of
- Maintains system sterility in both sides of the system via disconnect process
- Ensure simultaneous and secure removal of both membranes

Eliminate the need for tube sealers post disconnect process

Eliminates accidental disconnects

Enables broader application compatibility

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-GDC

Scan code to visit webpage



cpcworldwide.com/AseptiQuik-GDC



ASEPTIQUIK® G DC SERIES DIMENSIONS

POLYCARBONATE - For gamma irradiation applications.



TERMINATION	TUBING	METRIC EQ.	PART NO.	Δ	в	6
1/4" HOSE BARB Aseptiquik G DC Body	1/4"	6.4 mm	AQGDC17004	2.62" (66.5 mm)	4.37" (111.0 mm)	0.60" (15.2 mm)
3/8" HOSE BARB Aseptiquik g DC Body	3/8"	9.5 mm	AQGDC17006	2.62" (66.5 mm)	4.37" (111.0 mm)	0.60" (15.2 mm)
1/2" HOSE BARB Aseptiquik g DC Body	1/2"	12.7 mm	AQGDC17008	2.62" (66.5 mm)	4.66" (118.4 mm)	0.89" (22.6 mm)

PRODUCT DIMENSIONS









OPEN FORMAT CONNECTION TECHNOLOGY

WHAT IS OPEN FORMAT CONNECTION TECHNOLOGY?

CPC's open format connectors (also known as quick connectors or quick disconnectors) are a straightforward and simple way to incorporate media transfer technology between your single-use systems. Our open format products feature male and female connector halves with caps and plugs to seal off the fluid pathway of the single-use system. Simply removing the cap and plug from each half of the system and joining the male and female connection links the fluid pathways of your two separate systems.

CPC's open format products—the MPC Series, MPX Series, and MPU Series Connectors enable connections from 1/8" ID tubing to 1" ID tubing. In addition, the MPC and MPX connectors feature an ergonomic thumb latch and optional locking sleeve to prevent accidental disconnection.

Sanitary adapters within the MPC and MPX product lines facilitate integration of components into single-use or hybrid (single-use to stainless) process systems.

Back-to-back adapters within the MPC and MPX product lines can be used to connect two identical body components or insert components. For additional flexibility, reducer options enable connection between an MPC and MPX product to link 1/8" to 1/2" ID tubing.



FEATURES	BENEFITS
Ergonomic thumb latch \longrightarrow	Easy to operate – even with gloved hands
Parting line-free hose barb \longrightarrow	Creates seamless connection to tubing
Optional locking sleeve (MPC, MPX) \longrightarrow	Prevents accidental disconnection
Mix and match termination sizes	Enables flexibility in your application

NOMINAL FLOW PATH (NFP) SIZE



\bigcirc								
	NFP	CV RANGE						
МРХ	(0.5")	4.0 - 17						
BACK-TO-BACK Bodies & Inserts	(0.375" – 0.5")	23 - 32						



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the NFP.

OPEN FORMAT CONNECTION OPTIONS

	TERMINATIONS								PRODUCTFLOW-PATHFEATURESMATERIALS			STERILIZATION FLOW (CHOOSE 1) INFORMATION)W MATION	
	HOSE BARB				SANITARY		And		E. E.	ANN	FOW.	Complete Market	Helt Home	it when the second		
	1.81	llan	301	1/2	and an	"/	341	1.1.1. 	1000 1000 1000	SHIP	1000	N. OO	SPANNIE SPANNIE	A A A A A A A A A A A A A A A A A A A	ANN	C. M.
MPC	\checkmark	\checkmark	\checkmark						\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	1/4"	0.1 - 8.0
MPC SANITARY							\checkmark	\checkmark		\checkmark		\checkmark	\checkmark	\checkmark	1/4"	0.1 - 8.0
MPC BACK-TO-BACK Bodies											~	\checkmark	\checkmark	\checkmark	1/4" - 1/2"	2.0 - 6.0
MPC BACK-TO-BACK INSERTS										\checkmark		\checkmark	\checkmark	\checkmark	1/4" - 1/2"	2.0 - 6.0
MPX			\checkmark	\checkmark					\checkmark	\checkmark	~	\checkmark	\checkmark	\checkmark	1/2"	4.0 - 17
MPX SANITARY							\checkmark	\checkmark		\checkmark		\checkmark	\checkmark	\checkmark	1/2"	4.0 - 17
MPX BACK-TO-BACK Bodies											\checkmark	\checkmark	\checkmark	\checkmark	1/4" - 1/2"	23 - 32
MPX BACK-TO-BACK Inserts										\checkmark		\checkmark	\checkmark	\checkmark	1/4" - 1/2"	23 – 32
MPU					\checkmark	\checkmark				\checkmark		\checkmark	\checkmark	\checkmark	1"	18-41

MPC/MPX ASSEMBLY PROCEDURE





NOTE: If using a body or cap component with a locking sleeve, twist to the "locked" position after connecting the two halves to prevent accidental disconnection. When you are ready to disconnect once again, twist the locking sleeve to the "unlocked" position.

OPEN FORMAT CONNECTION TECHNOLOGY

MPC Series connectors

MPC Series Connecters add ease of use and security to critical fluid handling applications. Choose from a full line of connectors and configurations, including pressure sealing caps and plugs, in sizes to fit 1/8" to 3/8" tubing. MPC couplings offer optional locking sleeves to further guard against accidental disconnects. In addition, coupling halves can be rotated when connected to reduce tube kinks.

SPECIFICATIONS

OPERATING PRESSURE Vacuum to 60 psi, 4.1 bar

OPERATING TEMPERATURE

Polycarbonate: -40°F to 250°F (-40°C to 121°C) Polysulfone: -40°F to 300°F (-40°C to 149°C)

STERILIZATION

Gamma: Up to 50 kGy irradiation Autoclave: Polycarbonate: Up to 250°F (121°C),

30 minutes, up to 10 repetitions Sterilize uncoupled only **Polysulfone:** Up to 270°F (132°C), 60 minutes, up to 25 repetitions Sterilize uncoupled only

TERMINATIONS 1/8" to 3/8" ID (3.2mm to 9.5mm)

MATERIALS

Main components: Polycarbonate (purple tint) Polysulfone (amber tint) Locking sleeves: Polysulfone (white) Thumb Latches: Polycarbonate (white) PVDF (white) O-rings: Silicone (clear), platinum-cured

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES

Ergonomic thumb latch

Parting line-free hose barb –

Optional locking sleeve

TYPICAL FLOW RATE:

Cv Value Range: 0.1 - 8 for MPC hose barb terminations



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

NOTE

BENEFITS

Validation and Extractables data can be requested at cpcworldwide.com/MPC

DID YOU KNOW -

Easy to operate – even with gloved hands

Prevent potential leak path

Prevents accidental disconnection

Better flexibility to fit more applications

The MPC and MPX connectors are perfect for smaller bag systems for aliquoted media or other product stored in bags.

Scan code to visit webpage



cpcworldwide.com/MPC



MPC SERIES DIMENSIONS

COUPLING BODIES



5	TERMINATION	METRIC EQ.	POLYCARBONATE Part no.	POLYSULFONE Part no.	٥	B
2	1/8" HOSE BARB	3.2 mm ID	MPC17002T03	MPC17002T39	0.96" (24.4 mm)	1.10" (27.9 mm)
	1/4" HOSE BARB	6.4 mm ID	MPC17004T03	MPC17004T39	0.96" (24.4 mm)	1.30" (33.0 mm)
P	3/8" HOSE BARB	9.5 mm ID	MPC17006T03	MPC17006T39	0.96" (24.4 mm)	1.30" (33.0 mm)
	1/8" HOSE BARB W/ LOCK	3.2 mm ID	MPCK17002T03	MPCK17002T39	1.02" (25.9 mm)	1.10" (27.9 mm)
	1/4" HOSE BARB W/ LOCK	6.4 mm ID	MPCK17004T03	MPCK17004T39	1.02" (25.9 mm)	1.30" (33.0 mm)
	3/8" HOSE BARB W/ LOCK	9.5 mm ID	MPCK17006T03	MPCK17006T39	1.02" (25.9 mm)	1.30" (33.0 mm)

COUPLING INSERTS

TERMINATION	METRIC EQ.	POLYCARBONATE Part no.	POLYSULFONE Part no.	A	B
1/8" HOSE BARB	3.2 mm ID	MPC22002T03M	MPC22002T39M	0.60" (15.2 mm)	1.09" (27.7mm)
1/4" HOSE BARB	6.4 mm ID	MPC22004T03M	MPC22004T39M	0.60" (15.2 mm)	1.30" (33.0 mm)
3/8" HOSE BARB	9.5 mm ID	MPC22006T03M	MPC22006T39M	0.60" (15.2 mm)	1.30" (33.0 mm)

2

SEALING	COMP	ONENTS
----------------	------	---------------



SEALING CAP	W/LOCK	MATERIAL	A	B
MPC32003	MPCK32003	Polycarbonate	0.96" (24.4 mm)	1.30" (33.0 mm)
MPC32039	MPCK32039	Polysulfone	0.99" (25.2 mm)	1.30" (33.0 mm)
SEALING PLUG		MATERIAL	A	ß
SEALING PLUG MPC30003M		MATERIAL Polycarbonate	A 0.75" (19.1 mm)	B 1.24" (31.5 mm)

PRODUCT DIMENSIONS





B = Total Length

OPEN FORMAT CONNECTION TECHNOLOGY

MPX SERIES CONNECTORS

MPX Series Connectors add ease of use and security to your media transfer applications. Choose from a full line of connectors and configurations, including pressure sealing caps and plugs in sizes to fit 3/8" and 1/2" tubing. MPX couplings offer optional locking sleeves to further guard against accidental disconnects. In addition, coupling halves can be rotated when connected reducing tube kinks.

SPECIFICATIONS

OPERATING PRESSURE Vacuum to 60 psi, 4.1 bar

OPERATING TEMPERATURE

Polycarbonate: -40°F to 250°F (-40°C to 121°C) Polysulfone: -40°F to 300°F (-40°C to 149°C)

STERILIZATION

Gamma: Up to 50 kGy irradiation Autoclave: Polycarbonate: Up to 250°F (121°C),

30 minutes, up to 10 repetitions Sterilize uncoupled only **Polysulfone:** Up to 270°F (132°C), 60 minutes, up to 25 repetitions Sterilize uncoupled only

TERMINATIONS 3/8" to 1/2" ID (9.5mm to 12.7mm)

MATERIALS

Main components: Polycarbonate (purple tint) Polysulfone (amber tint) Locking sleeves: PVDF (white) Thumb Latches: Polycarbonate (white) PVDF (white) O-rings:

Silicone (clear), platinum-cured

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES

Ergonomic thumb latch

Parting line-free hose barb —

Optional locking sleeve –

Mix and match termination sizes —



BENEFITS

\rightarrow	Easy to operate – even with gloved hands
\rightarrow	Prevents potential leak path
\rightarrow	Prevents accidental disconnection
\rightarrow	Enables flexiblity in your system and/or application

TYPICAL FLOW RATE:

Cv Value Range: 4 - 17 for MPX



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/MPX

DID YOU KNOW

The MPC and MPX connectors are perfect for smaller bag systems for aliquoted media or other product stored in bags.

Scan code to visit webpage



cpcworldwide.com/MPX

MPX SERIES DIMENSIONS

COUPLING BODIES



TERMINATION	METRIC EQ.	FLOW	POLYCARBONATE Part no.	POLYSULFONE Part no.	۵	8
1/2" ID HOSE BARB	12.7 mm ID	.50"	MPX17803	MPX17839	1.28" (32.5 mm)	1.96" (49.8 mm)
1/2" ID HOSE BARB W/ LOCK	12.7 mm ID	.50"	MPXK17803	MPXK17839	1.28" (32.5 mm)	1.96" (49.8 mm)

COUPLING INSERTS



TERMINATION	METRIC EQ.	FLOW	POLYCARBONATE Part no.	POLYSULFONE Part no.	•	B
3/8" HOSE BARB	9.5 mm ID	.38"	MPX22603M	MPX22639M	0.85" (21.6 mm)	1.90" (48.3 mm)
1/2" HOSE BARB	12.7 mm ID	.50"	MPX22803M	MPX22839M	0.85" (21.6 mm)	1.90" (48.3 mm)

SEALING COMPONENTS



SEALING CAP	W/LOCK	MATERIAL	A	B
MPX32003	MPXK32003	Polycarbonate	1.28" (32.5 mm)	1.67" (42.4 mm)
MPX32039	MPXK32039	Polysulfone	1.28" (32.5 mm)	1.67" (42.4 mm)
SEALING PLUG		MATERIAL	Δ	B
MPX30003M		Polycarbonate	1.10" (27.9 mm)	1.66" (42.2 mm)
MPX30039M		Polysulfone	1.10" (27.9 mm)	1.66" (42.2 mm)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

PRODUCT DIMENSIONS



OPEN FORMAT CONNECTION TECHNOLOGY

MPC/MPX BACK-TO-BACK Series Adapters

MPC/MPX Back-to-Back Adapters give end users the flexibility of connecting single-use systems that feature identical coupling connections at the end of their tubing. Combining both MPC and MPX couplings provides a reducing option for users who need to transition between tubing diameters ranging from 1/8" to 1/2".

SPECIFICATIONS

OPERATING PRESSURE Vacuum to 60 psi, 4.1 bar

OPERATING TEMPERATURE

Polycarbonate: -40°F to 250°F (-40°C to 121°C) Polysulfone: -40°F to 300°F (-40°C to 149°C)

STERILIZATION

Gamma: Up to 50 kGy irradiation Autoclave: Polycarbonate: Up to 250°F (121°C), 30 minutes, up to 10 repetitions Sterlize uncoupled only Polysulfone: Up to 270°F (132°C), 60 minutes, up to 25 repetitions Sterlize uncoupled only

MATERIALS

Main Components:

Polycarbonate (purple tint) Polysulfone (amber tint)

Thumb Latches: Polycarbonate (white) PVDF (white)

O-rings: Silicone (clear), platinum-cured

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES

Compatible with MPC	
and MPX Series inserts	

Tubing reduction option

Ergonomic thumb latches

BENEFITS

Easy conversion to industry standard connections or single-use systems

Allows easy transition between multiple size tubing from 1/8" to $1/2"~{\rm ID}$

Easy to operate - even with gloved hands

TYPICAL FLOW RATE:

Cv Value Range: 2 - 32 for Back-to-Back hose barb terminations



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/Back-to-Back-Adapters

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/Back-to-Back-Adapters

DID YOU KNOW

The MPC and MPX connectors are perfect for smaller bag systems for aliquoted media or other product stored in bags.



MPC/MPX BACK-TO-BACK SERIES DIMENSIONS

MPC/MPX BACK-TO-BACK INSERT ADAPTERS - Polysulfone

	PART NO.	ТҮРЕ	A	B
	MPC22C2239M	MPC to MPC	0.74" (18.8)	2.04" (51.0)
100 M	MPC22X2239M	MPC to MPX	0.98" (25.0)	2.42" (61.5)
200	MPX22X2239M	MPX to MPX	0.98" (25.0)	2.73" (69.5)

MPC/MPX BACK-TO-BACK BODY ADAPTERS - Polycarbonate

	PART NO.	ТҮРЕ	A	8
The start	MPC17C1703	MPC to MPC	0.96" (24.5)	1.81" (46.0)
A Share	MPX17X1703	MPX to MPX	1.28" (32.5)	2.44" (62.0)
Character	MPC17X1703	MPC to MPX	1.28" (32.5)	2.13" (54.1)

MPC/MPX BACK-TO-BACK BODY ADAPTERS - Polysulfone

	PART NO.	ТҮРЕ	۵	B	
A M	MPC17C1739	MPC to MPC	0.96" (24.5)	1.81" (46.0)	
A	MPX17X1739	MPX to MPX	1.28" (32.5)	2.44" (62.0)	
Cha	MPC17X1739	MPC to MPX	1.28" (32.5)	2.13" (54.1)	

PRODUCT DIMENSIONS







NOTES -

OPEN FORMAT CONNECTION TECHNOLOGY

MPC/MPX SANITARY Series connectors

MPC/MPX Sanitary Connectors attach directly to 3/4", 1" and 1-1/2" sanitary terminations to provide greater flexibility for integrating components into single-use or hybrid (single-use to stainless) process systems. Standard bag systems with quick couplings can be easily connected to equipment with sanitary terminations, while single-use cartridge filters can be converted to incorporate quick couplings for greater system modularity.



SPECIFICATIONS

OPERATING PRESSURE: Vacuum to 60 psi, 4.1 bar

OPERATING TEMPERATURE: -40°F to 300°F (-40°C to 149°C)

STERILIZATION:

Gamma: Up to 50 kGy irradiation **Autoclave:** Up to 270°F (132°C) for 60 minutes, up to 25 repetitions. Sterilize uncoupled only.

TERMINATIONS: 3/4", 1" and 1-1/2" sanitary

MATERIALS:

Main components: Polysulfone (amber tint) Thumb Latches: PVDF (white) O-rings: Silicone (clear), platinum-cured

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

FEATURES

3/4", 1" and 1-1/2	" sanitary terminations	
--------------------	-------------------------	--

Compatible with MPC	
and MPX Series couplings	

Integral coupling adapters -

ADCF-free materials

TYPICAL FLOW RATE:

Cv Value Range: 3 - 17



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/Sanitary

BENEFITS

 \rightarrow

 \rightarrow

Install to equipment with sanitary gaskets and sanitary clamps

- Quick and easy connections to industry standard plastic couplings
- Provides flexibility to easily convert sanitary terminations on filter cartridge or equipment

Meet BSE/TSE requirements

NOTE

Validation and Extractables data can be requested on at cpcworldwide.com/Sanitary

DID YOU KNOW

MPC and MPX Sanitary connectors provide greater flexibility for filter installation.

MPC/MPX SANITARY SERIES DIMENSIONS

COUPLING BODIES - Polysulfone



PART NO.	SIZE	A	B	С
MPC3301239	3/4"	0.98" (24.9)	1.40" (35.6)	1.0" (25.4)
MPX3301239	3/4"	1.28" (32.5)	1.70" (43.2)	1.0" (25.4)
MPC3301639	1"	1.50" (38.1)	1.40" (35.6)	1.50" (38.1)

COUPLING INSERTS - Polysulfone



PART NO.	SIZE	A	B
MPC44012T39M	3/4"	0.98" (24.9)	1.40" (35.6)
MPC44024T39M	1-1/2"	1.98" (50.3)	1.40" (35.6)
MPX44012T39M	3/4"	0.98" (24.9)	1.71" (43.4)
MPX44024T39M	1-1/2"	1.98" (50.3)	1.71" (43.4)

PRODUCT DIMENSIONS







NOTES -

OPEN FORMAT CONNECTION TECHNOLOGY

MPU Series connectors

MPU Connectors' twist-to-connect design

features an easy-to-use locking mechanism that guards against accidental disconnects and provide a reliable, secure connection. The 3/4" and 1" hose barbs provide smooth, rapid media transfer.



SPECIFICATIONS

OPERATING PRESSURE Vacuum to 35 psi, 2.4 bar

OPERATING TEMPERATURE -40°F to 300°F (-40°C to 149°C)

STERILIZATION

Gamma: Up to 50 kGy irradiation Autoclave: Up to 270°F (132°C), 60 minutes, up to 25 repetitions Sterilize uncoupled only

TUBING SIZE 3/4" ID (19.0 mm), 1" ID (25.4 mm)

MATERIALS

Main components: Polysulfone (amber tint) O-rings:

Silicone (clear), platinum-cured

WARNING: Pressure, temperature, chemicals, and operating environment can affect the performance of couplings. It is the customer's responsibility to test the suitability of CPC's products in their own application conditions.

Scan code to visit webpage



cpcworldwide.com/MPU

FEATURES

3/4" and 1" hose barb

Locking feature

Sharp barb end

Shrouded, leak-free seal & smooth, _____ internal flow path

TYPICAL FLOW RATE:

Cv Value Range: 18 - 41 for MPU



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

BENEFITS

Facilitates rapid fill and empty of bioprocessing bags

Guards against accidental disconnects

Minimizes fluid turbulence and dead space

Protect valuable fluids and eliminate potential to contaminate fluid path

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/MPU

DID YOU KNOW

The MPU connectors are perfect for attaching to single-use mixers or single-use bioreactors when a large amount of media needs to be transferred.

MPU SERIES DIMENSIONS

COUPLING BODIES - Polysulfone

TERMINATION	METRIC EQ.	PART NO.	A	в
3/4" ID HOSE BARB	19.0 mm ID	MPU171239	1.75" (44.5 mm)	2.37" (60.2 mm)
1" ID HOSE BARB	25.4 mm ID	MPU171639	1.75" (44.5 mm)	2.37" (60.2 mm)

COUPLING INSERTS - Polysulfone

0	TERMINATION 3/4" ID HOSE BARB	METRIC EQ. 19.0 mm ID	PART NO. Mpu221239M	▲ 1.9" (48.3 mm)	© 2.88" (73.2 mm)
	1" ID HOSE BARB	25.4 mm ID	MPU221639M	1.9" (48.3 mm)	2.88" (73.2 mm)

SEALING COMPONENTS - Polysulfone

SEALING CAP	MATERIAL	A	8
MPU32039	Polysulfone	1.75" (44.5mm)	.79" (20.1mm)
MPU30039M	Polysulfone	1.56" (39.6mm)	1.38" (35.1mm)

PRODUCT DIMENSIONS



MPU ASSEMBLY PROCEDURE





STEAM-IN-PLACE (SIP) TECHNOLOGY

WHAT IS STEAM-IN-PLACE (SIP) TECHNOLOGY?

The innovative three port design of CPC's Steam-Thru® technology enables uncomplicated actuation between the SIP steam pathway and fluid transfer (flow) pathway. This provides an easy and sterile connection between flexible single-use tubing and stainless-steel processing equipment. The Steam-Thru product lines offer versatility within your process by having connection options for 3/8" and 1/2" ID tubing on the single-use side of the connector and 3/4" and 1-1/2" sanitary connections for attachment to the stainless-steel processing equipment.



FEATURES

BENEFITS

Innovative three-port design \longrightarrow	Allows for a true steam-through SIP process to eliminate dead legs and the need for laminar flow hoods
Intuitive valve design \longrightarrow	Enables sterile connection/disconnection while permitting a high media flow rate
Steam-Thru II Thumb Latch \longrightarrow	Secures valve position, provides visual indicator of process stage
3/4" and 1-1/2" Sanitary Terminations \longrightarrow	Easily connects to process equipment
Additional features and benefits for the Asep	tiQuik Steam-Thru Combination:
Genderless AseptiQuik \longrightarrow	Eases single-use systems specifications with one part number for both halves
FLIP-CLICK-PULL	Intuitive three-step connection process reduces risk of operator error
$\operatorname{CPC}\operatorname{Click} \longrightarrow$	Audible confirmation of assembly with no additional hardware required
Sanitary interface between the two connectors \longrightarrow	More secure connection than tubing with cable ties

NOMINAL FLOW PATH (NFP) SIZE



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the NFP.

CPC

STEAM-THRU CONNECTION OPTIONS

	; Te	SINGL Sys Rmin	E USE Tem Iatioi	E NS	EQUIF Po Termii	PMENT RT Nation	CONDENSATE Port Termination	PR Fe	ODU Atur	ICT Res	STERILI Opti (Cho(ZATION ONS DSE 1)	SIP	FLO Inform	FLOW INFORMATION	
	HO Ba	SE RB	ASE Qu	PTI- Iik	SANI	TARY	SANITARY	 				All			HI HI	
	3,0,	1102	405	400	8 1 1 1 1		100 m	Sterner -	HI HI HI HI		the state of the s			Annual Contraction of the second seco		
STEAM-THRU	\checkmark	\checkmark			\checkmark	\checkmark	\checkmark	~			\checkmark	\checkmark	\checkmark	1/2"	4.0 - 7.0	
STEAM-THRU II	\checkmark	\checkmark			\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	1/2"	3.0 - 9.0	
ASEPTIQUIK S Steam-thru II			\checkmark		\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	3/8"	1.0 - 2.0	
ASEPTIQUIK G Steam-thru II				\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	1/2"	3.0 - 9.0	

STEAM-THRU II ASSEMBLY PROCEDURE



CPC

STEAM-THRU® SERIES CONNECTORS

Steam-Thru® Connectors allow a quick and easy sterile connection between stainless steel biopharmaceutical processing equipment and disposable bag and tube assemblies. The single-use design saves time and money by eliminating unnecessary cleaning procedures and reducing validation burden associated with reusable components.

SPECIFICATIONS

OPERATING CONDITIONS (Fluid Transfer)

STEAM POSITION

Temperature: 39°F to 104°F (4°C to 40°C) Pressure: Up to 30 psi, 2.1 bar (Steam-Thru) Up to 35 psi, 2.4 bar (Steam-Thru II)

FLOW POSITION

Temperature: 39°F to 104°F (4°C to 40°C) **Pressure:** Vacuum to 20 psi, 1.4 bar

STERILIZATION

Full Connector Assembly Gamma: Up to 50 kGy irradiation Autoclave: Up to 265°F (129°C) for 60 minutes

Steam-In-Place (SIP): Up to 266°F/130°C (up to ~24psi), 60 minutes (Steam-Thru) Up to 275°F/135°C (up to ~31psi), 60 minutes (Steam-Thru II)

TERMINATIONS

3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb (Steam-Thru) 3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb and 3/4" sanitary (Steam-Thru II)

MATERIALS

Connection: Polysulfone (amber tint) **O-rings:** Silicone (clear), platinum-cured **Removeable Sleeve:** Polycarbonate

NOTE: Steam pressures are estimated based upon information in Steam Tables found in literature

FEATURES

Innovative three-port design _____

Valve design _____

Steam-Thru II thumb latch -

3/4" and 1-1/2" sanitary terminations -

TYPICAL FLOW RATE:

Cv Value Range: 4.0 - 7.0 for Steam-Thru 3.0 - 9.0 for Steam-Thru II



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

Scan code to visit webpage



cpcworldwide.com/STC

BENEFITS

Allows a true steam-through SIP process which eliminates "dead legs" and the need for laminar flow hoods

Allows sterile connection and disconnection and permits high media flow rate

Secures valve position, provides visual indicator of process stage

Easily connects to process equipment

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/STC

DID YOU KNOW

Steam-Thru connectors are perfect for any hybrid processing at your facility. If using a stainless-steel bioreactor and single-use systems, easily make a sterile connection between the two systems with a Steam-Thru connector on your single-use system and mounting it directly onto your bioreactor.

STEAM-THRU SERIES DIMENSIONS

POLYSULFONE with polycarbonate sleeve



TERMINATION PART NO. A B G D e 3/4" X 3/4" SANITARY X 1/2" HB STC1700000 1.20" (30.5 mm) 5.09" (129.3 mm) 4.44" (112.8 mm) 2.00" (50.8 mm) 0.89" (22.6 mm) 3/4" X 3/4" SANITARY X 3/8" HB STC1700100 0.60" (15.2 mm) 1.20" (30.5 mm) 4.80" (121.9 mm) 4.15" (105.4 mm) 2.00" (50.8 mm) 3/4" X 1-1/2" SANITARY X 1/2" HB STC1700200 1.20" (30.5 mm) 5.09" (129.3 mm) 4.44" (112.8 mm) 2.00" (50.8 mm) 0.89" (22.6 mm) 3/4" X 1-1/2" SANITARY X 3/8" HB STC1700300 1.20" (30.5 mm) 4.80" (121.9 mm) 4.15" (105.4 mm) 2.00" (50.8 mm) 0.60" (15.2 mm)

STEAM-THRU CONFIGURATIONS

Steam-Thru Connection's three-port design allows steam to pass directly through the lower ports to "steam on" to stainless equipment. After the SIP cycle is completed, the connector's valve is actuated, creating a sterile flow path to single-use systems.



STEAM-THRU II SERIES DIMENSIONS

POLYSULFONE with polycarbonate sleeve



TERMINATION	PART NO.	A	B	0	D	•
3/4" X 3/4" SANITARY X 1/2" HB	STC2020000	1.42" (36.1 mm)	6.84" (173.7 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.88" (22.4 mm)
3/4" X 3/4" SANITARY X 3/8" HB	STC2020100	1.42" (36.1 mm)	6.76" (171.7 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.80" (20.3 mm)
3/4" X 1-1/2" SANITARY X 1/2" HB	STC2020200	1.42" (36.1 mm)	6.84" (173.7 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.88" (22.4 mm)
3/4" X 1-1/2" SANITARY X 3/8" HB	STC2020300	1.42" (36.1 mm)	6.76" (171.7 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.80" (20.3 mm)
3/4" X 3/4" SANITARY X 3/4" Sanitary	STC2020900	1.42" (36.1 mm)	6.60" (167.6 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.62" (15.7 mm)
3/4" X 1-1/2" SANITARY X 3/4" Sanitary	STC2021000	1.42" (36.1 mm)	6.60" (167.6 mm)	5.93" (150.6 mm)	2.40" (61.0 mm)	0.62" (15.7 mm)

All measurements are in inches (millimeters) unless otherwise noted. Tubing must meet stated inside and outside diameters.

STEAM-THRU II CONFIGURATIONS

Steam-Thru II Connections offer the flexibility of "steam on" and "steam off" functionality. The innovative design allows the valve to be returned to the steam position enabling a second SIP cycle following media transfer. The "steam off" disconnection of single-use systems minimizes cross-contamination risks associated with reusable components.



ASEPTIQUIK® STC II SERIES CONNECTORS

AseptiQuik® STC Connectors combine the AseptiQuik® sterile connector and the Steam-Thru® II SIP connector, giving manufacturers greater flexibility between hybrid stainless steel and single-use processing equipment.

The union of the two connectors into a single unit through a sanitary clamp allows an AseptiQuik sterile connection to be steamed on to stainless equipment via SIP. After the SIP cycle, a wide range of single-use systems can be connected.



SPECIFICATIONS

OPERATING CONDITIONS

(Fluid Transfer) Steam Position: Temperature: Up to 275°F (135°C) for 60 minutes

PRESSURE

Up to 35 psi, 2.4 bar Flow Position: Temperature: 39°F to 104°F (4°C to 40°C) Pressure: Up to 20 psi, 1.4 bar

STERILIZATION

Full Connector Assembly: Gamma: Up to 50 kGy irradiation Autoclave: High Temp (HT) Version Up to 266°F (130°C) for 60 minutes (AQGSTC) Steam-In-Place (SIP):

Up to 275°F (135°C) (~31 psi) for 60 minutes

TERMINATIONS

3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb (Steam-Thru) 3/8" (9.5mm) to 1/2" (12.7mm) ID hose barb and 3/4" sanitary (Steam-Thru II)

MATERIALS:

Main Components: AseptiQuik - Polycarbonate (white) Steam-Thru II - Polysulfone (amber tint) Seals: Silicone (clear), platinum-cured Removable Sleeve: Polycarbonate (white) Pull Tabs:

Polycarbonate (blue, standard version) Polycarbonate (white, HT version)

Membrane:

Polyethylene (standard version) Hydrophobic polyethersulfone (HT versions), PTFE strip sticker **Clamp:** Nylon 66 (white)

NOTE: Steam pressures are estimated based upon information in Steam Tables found in literature

Scan code to visit webpage



cpcworldwide.com/AseptiQuik-STC

TYPICAL FLOW RATE:

Cv Value Range: 1 - 2 for AQSSTCII 3 - 9 for AQGSTCII



Cv values represent the approximate expected flow rate in gallons per minute of water at room temperature for a 1 PSI pressure drop. The flow is generally constrained by the smallest diameter, which in some cases will be the termination diameter and not the Nominal Flow Path.

NOTE

Validation and Extractables data can be requested at cpcworldwide.com/AseptiQuik-STC

DID YOU KNOW

The AQSSTC provides the same sterile hybrid technology as the Steam-Thru II, but in an even more compact form. The AQSSTC has an AseptiQuik mounted on the single-use port of the Steam-Thru Connector. Meaning fewer single-use systems mounted on your bioreactor during the SIP process. Connect the other end of the AQS at any point after the SIP process and before actuating to the flow position.

ASEPTIQUIK STC SERIES DIMENSIONS

COMBINATION PRODUCT - Polycarbonate with blue pull tabs and Steam-Thru II- For gamma irradiation applications.



TERMINATION	PART NO.	A	B	C
3/4" X 3/4" SANITARY STEAM-THRU II With Aseptiquik & Body	AQGSTC2330900	3.05" (77.5 mm)	8.26" (209.8 mm)	7.36" (186.9 mm)
3/4" X 1-1/2" SANITARY STEAM-THRU II With Aseptiquik & Body	AQGSTC2331000	3.05" (77.5 mm)	8.26" (209.8 mm)	7.36" (186.9 mm)
3/4" X 3/4" SANITARY STEAM-THRU II With Aseptiquik S Body	AQSSTC2330900	2.45" (62.2 mm)	8.20" (208.2 mm)	7.30" (185.4 mm)
3/4" X 1-1/2" SANITARY STEAM-THRU II WITH ASEPTIQUIK S BODY	AQSSTC2331000	2.45" (62.2 mm)	8.20" (208.2 mm)	7.30" (185.4 mm)

COMBINATION PRODUCT - Polycarbonate HT with white pull tabs and Steam-Thru II- For autoclave or gamma irradiation applications.



TERMINATION	PART NO.	٥	B	G
3/4" X 3/4" SANITARY STEAM-THRU II With Aseptiquik G Body	AQGSTC2330900HT	3.05" (77.5 mm)	8.26" (209.8 mm)	7.36" (186.9 mm)
3/4" X 1-1/2" SANITARY STEAM-THRU II With Aseptiquik G Body	AQGSTC2331000HT	3.05" (77.5 mm)	8.26" (209.8 mm)	7.36" (186.9 mm)
3/4" X 3/4" SANITARY STEAM-THRU II With Aseptiquik S Body	AQSSTC2330900HT	2.45" (62.2 mm)	8.20" (208.2 mm)	7.30" (185.4 mm)
3/4" X 1-1/2" SANITARY STEAM-THRU II With Aseptiquik S Body	AQSSTC2331000HT	2.45" (62.2 mm)	8.20" (208.2 mm)	7.30" (185.4 mm)

PRODUCT DIMENSIONS



1 72	1 11	87			
Ŀ.	ы.	ΛL		ĸ	F
	l - i	m I.	U	11	LJ

BENEFITS

Genderless AseptiQuik \longrightarrow	Eases single-use systems specifications with one- part number for both halves
FLIP-CLICK-PULL	Intuitive three-step connection process reduces risk of operator error
CPC Click \longrightarrow	Audible confirmation of assembly with no additional hardware required
Innovative three-port steam design \longrightarrow	Allows a true steam-through SIP process which eliminates "dead legs"
Steam valve design \longrightarrow	Allows sterile connection and disconnection to stainless equipment and permits a high media flow rate
Sanitary interface between \longrightarrow the two connectors	More secure connection than tubing with cable ties



PRODUCT INDEX

ASEPTIC CONNECTION

TECHNOLOGY	.10-23
MICROCNX [®] Series Connectors	12
ASEPTIQUIK [®] S Series Connectors	16
ASEPTIQUIK [®] G Series Connectors	18
ASEPTIQUIK [®] L Series Connectors .	20
ASEPTIQUIK [®] W Series Connectors	22

STERILE DISCONNECTION

TECHNOLOGY										.24	-29
HFC Disconnec	t	Se	rie	es (Co	nn	ес	to	rs		. 28

ASEPTIC COMBINATION

CONNE	CTO	RS									.3	0-	35
ASEPT	IQUI	K® G	DC	Se	erie	S	Со	nn	ec	toı	ſS.		34

OPEN FORMAT CONNECTION

TECHNOLOGY
MPC Series Connectors 40
MPX Series Connectors
MPC/MPX Back-to-Back Series Adapters 44
MPC/MPX Sanitary Series Connectors 46
MPU Series Connectors

STEAM-IN-PLACE CONNECTION

FECH	NOL	00	iY				-		-		-		.5	0-	5	7
STE/	AM-1	THR	U®	Se	eri	es	Сс	nr	nec	cto	rs				5	4
ASEI	PTIQ	UIK	(® S	Y												







Colder Products Company®

COLDER PRODUCTS COMPANY 2820 CLEVELAND AVE. N., Roseville, minnesota 55113 U.S.A.

 PHONE:
 +1 (651) 645-0091

 FAX:
 +1 (651) 645-5404

 TOLL FREE:
 +1 (800) 519-7633

COLDER PRODUCTS COMPANY GMBH Kurhessenstrasse 15 64546 Mörfelden-Walldorf Germany

PHONE: +49 6105 9743-003

DOVER (SHANGHAI) INDUSTRIAL CO., LTD. Room 1802, Building A, IBC, 391 Guiping Road, Xuhui District, Shanghai, China 200233

PHONE: +86 21 2411 2666 **TOLL FREE:** +86 400 990 1978

cpcworldwide.com/bio info@cpcworldwide.com

CPC WARRANTY STATEMENT: CPC (Colder Products Company) warrants its products against defects in workmanship and materials for a period of 12 months from the date of sale by CPC to its initial customer (regardless of any subsequent sale of the products). This warranty is void if the product is misued, altered, tampered with or is installed or used in a manner that is inconsistent with CPC's written recommendations, specifications and/or instructions, or fails to perform due to normal wear and tear. CPC does not warrant the suitability of the product for any particular application. Determining product application suitability is solely the customer's responsibility. CPC is not liable for special, indirect, incidental, consequential or other damages including, but not limited to, loss, damage, personal injury, or any other expense directly or indirectly arising from the use of or inability to use its products either separately or in combination with other products. ALL OTHER WARRANTIES EXPRESS OR IMPLIED, WHETHER ORAL, WRITTEN OR IN ANY OTHER FORM, INCLUDING BUT NOT LIMITED TO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE EXPRESSLY EXCLUDED. The sole and exclusive remedy under this warranty is limited, at the option of CPC, to replacement of the defective product or an account credit in the amount of the original selling price. All allegedly defective CPC products must be returned prepaid transportation to CPC, together with information describing the product's application and performance, unless otherwise authorized in writing by CPC.

CPC PATENT STATEMENT: CPC takes pride in its innovative quick disconnect coupling and fittings solutions, many of which have been awarded United States and international patents. CPC has a strong tradition of leadership in the quick disconnect market, and aggressively pursues and protects its proprietary information and intellectual property. In cases where it is practical and has a benefit to its customers, CPC has licensed its proprietary technology. Please contact CPC to discuss your unique needs.

CPC TRADEMARK STATEMENT: AseptiQuik[®], BottleQuik[®], BreakAway[®], ChemQuik[®], DrumQuik[®], FitQuik[®], IdentiQuik[®], Nu-Seal[®], SnapQuik[®], Steam-Thru[®], Softube[®] are registered trademarks with the U.S. Patent & Trademark Office. All other trademarks or service marks are property of their respective owners. WARNING: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC. It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer.

COPYRIGHT © 2023 COLDER PRODUCTS COMPANY. All rights reserved. Colder Products Company,Colder Products and CPC are registered trademarks with the US Patent & Trademark Office.

WARNING: Due to the wide variety of possible fluid media and operating conditions, unintended consequences may result from the use of this product, all of which are beyond the control of CPC It is the user's responsibility to carefully determine and test for compatibility for use with their application. All such risks shall be assumed by the buyer.

