



Editorial Contact:
Patty Hoffman
Padilla Speer Beardsley
Tel: 612-455-1781

Email: phoffman@padillaspeer.com



For a high resolution image, visit:
www.colder.com/AseptiQuikS-X_photos

Colder Products Company to Showcase New Small and Large Single-Use Sterile Connectors at INTERPHEX 2012

New connectors provide biopharmaceutical manufacturers repeatable, reliable performance for greater process flexibility

St. Paul, Minn. – May 1, 2012 – [Colder Products Company](http://www.colder.com), the leader in single-use connection technology, is unveiling **AseptiQuik® S** small-format connectors and **AseptiQuik® X** large-format connectors at INTERPHEX 2012 (Booth #2537) at the Javits Center in New York. Biopharmaceutical manufacturers now have access to a full line of robust, sterile connectors with flow configurations ranging from 1/8 inch to 1-inch. Colder experts will be available to show event attendees how they can easily incorporate the new connection technologies into their single-use applications.

“Manufacturers will benefit from the wide range of sterile connection whether taking a process sample or harvesting a 1,000 liter mixing system,” said John Boehm, business unit manager, bioprocessing markets. “The new expanded sterile connector sizes demonstrate our continued commitment to developing solutions for greater process flexibility and improved efficiencies.”

The **AseptiQuik® S** small format 1/8 inch and 1/4 inch connectors offer connections for small-flow applications and feature a genderless, easy-to-use design for reliable performance without the need for clamps, fixtures or tube welders. Applications include single-use bioreactors, cell banking, and a range of sampling applications. Terminations will include 1/8 inch hose barb, 1/4 inch hose barb, 3/4 inch sanitary and MPC male insert.

The **AseptiQuik® X** large format 1-inch connector provides quick and easy sterile connections for high-flow applications, even in non-sterile environments. With this new connector, time to transfer large volumes of media is greatly reduced saving valuable time and improving process efficiencies.

(more)

Applications include large format single-use bioreactors, mixers and transfer lines. Terminations will include 3/4 inch hose barb, 1 inch hose barb and 1-1/2 inch sanitary.

Colder's bioprocessing connections are manufactured in an ISO Class 7 certified cleanroom and meet USP Class VI standards.

To learn more about how to benefit from the new AseptiQuik connectors or other Colder bioprocessing solutions, visit www.colder.com/bio. For information on any of the other 7,500 Colder connection solutions, please visit www.colder.com or call 800-444-2474.

#

About Colder Products Company

Colder Products Company is the leading provider of quick disconnect couplings and sterile connection technologies for the biopharmaceutical, medical and special industrial markets. For a free catalog or more information, contact: Colder Products Company, 1001 Westgate Drive, St. Paul, MN 55114. Phone: 651-645-0091. Toll-free: 800-444-2474. Fax: 651-645-5404. Web: www.colder.com.

CPC, Colder Products and Colder Products Company are registered trademarks with the U.S. Patent & Trademark office.

(more)