



Editorial Contact:

Patty Hoffman
Padilla Speer Beardsley
Tel: (612) 455-1781
Email: phoffman@padillaspeer.com



For a high-resolution image, visit:

www.colder.com/snapquikphoto

Colder Products Company Launches SnapQuik™ Connectors

Connectors Aim to Reduce Tubing Misconnections, Improve Safety in Medical Applications

St. Paul, Minn. – Feb. 14, 2012 – In medical care situations, there are as many as 40 device connectors used with a single patient. Despite the best efforts of well-trained, highly qualified healthcare professionals, medical-tubing misconnections continue to occur due to the widespread use of standard luer connectors. [Colder Products Company](#), a leading developer of secure, leak-free connectors, has developed the new [SnapQuik™](#) connectors to give medical personnel a safer choice for small-bore connections in air-driven devices and other medical applications.

The SnapQuik connectors are a cost-effective and significantly more reliable alternative to luer fittings. Colder's SnapQuik connectors do not mate with luer fittings, helping users avoid misconnections that can lead to unintentional patient injury or death.

“The compact and intuitive SnapQuik product line is designed to help reduce potentially deadly medical-tubing misconnections,” said Jim Brown, medical business unit manager for Colder Products Company. “These innovative connectors are also designed with a streamlined exterior and a unique internal-latching mechanism. This easy-to-clean design reduces the risk of dirt build-up and minimizes the potential for patient infection.”

SnapQuik connectors are ideal for applications such as air-inflated cuffs, surgical devices and monitoring equipment. The connectors deliver high-flow fluid transfer in a small profile. SnapQuik connectors are designed for 3/32” and 1/8” tubing sizes.

The Colder SnapQuik products deliver an audible “click” to assure reliable connections are made and feature breakaway latch mechanisms for quick and easy disconnections. The SnapQuik series is the latest offering in Colder’s platform of products for leak-free fluid connections.

Colder will unveil the SnapQuik product line at MD&M West (booth #2455), Feb. 14-16 at the Anaheim Convention Center in California. Attendees will have the opportunity to meet with Colder experts and discuss this new small-bore connector for applications requiring a design that eliminates the potential for misconnections associated with luer fittings.

#

About Colder Products Company

Colder Products Company is the leading provider of quick disconnect couplings, fittings and connectors to life sciences, industrial and chemical handling markets. For a free catalog or more information, contact: Colder Products Company, 1001 Westgate Drive, St. Paul, MN 55114. Phone: 651-645-0091. Toll-free: 800-444-2474. Fax: 651-645-5404. Web: www.colder.com.

CPC, Colder Products and Colder Products Company are registered trademarks with the U.S. Patent & Trademark office.